

# BioDEG

Biodegradation and Corrosion Simulation using Finite Element

User Manual

Version 0.8

(generated December 17, 2021)

Mojtaba Barzegari  
Liesbet Geris

[BioDeg website](#)

**KU LEUVEN**

---

Copyright (c) 2019-2021 University of Leuven and [BioDeg authors](#).

# Contents

<b>1</b>	<b>Introduction</b>	<b>3</b>
1.1	Authors . . . . .	3
1.2	Acknowledgments . . . . .	3
1.3	Referencing BIODEG . . . . .	3
<b>2</b>	<b>Useful background information</b>	<b>4</b>
<b>3</b>	<b>Installation</b>	<b>4</b>
3.1	Easy installation . . . . .	4
3.2	Advanced installation for improved performance/flexibility . . . . .	5
3.2.1	Compiling and installing external libraries . . . . .	5
3.2.1.1	PETSc and Qt . . . . .	5
3.2.1.2	FreeFEM . . . . .	5
3.2.2	Building and installing BIODEG . . . . .	7
3.2.2.1	Build BIODEG UI using Qt Creator IDE . . . . .	7
3.2.2.2	Build BIODEG UI using Qt tools . . . . .	7
<b>4</b>	<b>Running BioDeg</b>	<b>8</b>
4.1	Configuring the simulation . . . . .	8
4.1.1	Example 1 . . . . .	8
4.1.2	Example 2 . . . . .	8
<b>5</b>	<b>Future extensions to BioDeg</b>	<b>8</b>
<b>6</b>	<b>Finding answers to more questions</b>	<b>8</b>
<b>A</b>	<b>Run-time input parameters</b>	<b>8</b>
A.1	Geometry and mesh parameters . . . . .	9
A.2	Materials and boundary conditions parameters . . . . .	11
A.3	Solver parameters . . . . .	14
A.4	Output parameters . . . . .	16
	<b>Index of run-time parameters with section names</b>	<b>19</b>

# 1 Introduction

BIODEG is an open source software written in FreeFEM (a domain-specific language for finite element programming), C++, and Python for modeling the degradation of metallic biomaterials and simulating the biodegradation behavior of medical devices, implants and scaffolds in corrosion experiments. It can handle any geometry of desire and supports parallel computing to simulate large scale models.

## 1.1 Authors

BIODEG is developed by the [Biomechanics Research group at KU Leuven and University of Liege](#). The code is currently maintained by its principal developer, who manage the development of the mathematical models and the core functionalities.

### Principal developer

- Mojtaba Barzegari (University of Leuven, Belgium)

### Previous Contributors

- Yann Guyot (University of Liege, Belgium)
- Piotr Bajger (University of Oxford, UK)

### Mentor

- Liesbet Geris (University of Leuven, Belgium)

### Chemist contributors

(who has helped to validate the models)

- Sviatlana V. Lamaka (Helmholtz-Zentrum Hereon, Germany)
- Di Mei (Zhengzhou University, China)
- Cheng Wang (Helmholtz-Zentrum Hereon, Germany)

## 1.2 Acknowledgments

The development of BIODEG open source code is financially supported by the Prosperos project, funded by the Interreg VA Flanders – The Netherlands program, CCI grant no. 2014TC16RFCB046 and by the Fund for Scientific Research Flanders (FWO), grant G085018N. The developers also acknowledge support from the European Research Council under the European Union's Horizon 2020 research and innovation programmes, ERC CoG 772418.

## 1.3 Referencing BioDeg

Please refer to [BIODEG repository](#), section "Publications and referencing" to properly cite the use of BIODEG in your scientific work.

## 2 Useful background information

You may refer to the following articles for a background of the methods and algorithms implemented in BIODEG.

In addition, below are some useful references on finite element method and some online resources that provide a background of finite elements and their application to the solution of partial differential equations.

## 3 Installation

Installing BIODEG is a straightforward procedure. You need to install a couple of prerequisites and download and run BIODEG; that's all you need to do. But for advanced users, it might be necessary or more interesting to build everything from scratch to have more control on customizing features and improving performance. As a result, we have provided 2 sets of instructions, one for easy installation using the compiled binaries and one for building things from the source codes. The installation instructions are provided for Linux and Windows operating systems, but the procedure should be very similar for macOS.

It is possible to use BIODEG without the user interface (UI) if this scenario is required by the user (like for running it on a super-computer). The core of BIODEG is written in FreeFEM, so in this case, all you need to do is to installing/building FreeFEM and the required libraries, and then, cloning the [BIODEG core repository](#) and running the code according the provided instruction in the README file of the repository.

The BIODEG UI contains all the bundles for pre-processing and post-processing simulation input/results. These features are being hosted on their own repositories (), but with obtaining the user interface, you can have them all together. If you choose to use BIODEG without the user interface and still want to use the provided script for pre/post-processing, you may need to obtain them separately.

For building the source codes, we assume standard tools and libraries like CMake, compilers (for C, C++ and Fortran), and MPI libraries are pre-installed on your machine. If you are going to build BIODEG and required dependencies on a super-computer or cluster, you should notice that most high-performance computers would have the latest version of these compilers and libraries in the default environment.

### 3.1 Easy installation

The simplest way to install and run BIODEG is via the pre-built binaries you can download from GitHub. The same principle applies to prerequisites, which is in this case FreeFEM only. So, following these steps will install BIODEG on your machine:

1. Download FreeFEM installer for the platform you use and install it. You can find the .exe installer for Windows and the .deb installer for Linux (Ubuntu) in the [Release page of FreeFEM repository](#). You will find these files under the Assets section of the latest (or any other version). Execute the download file and follow the installation procedure appearing on your screen.
2. Download BIODEG tarballs for your preferred platform (Windows or Linux) from the [Release page of BIODEG repository](#). This is indeed the BIODEG UI bundle that contains the BIODEG core, the user

interface, the pre-processor, and the post-processor.

3. Extract the downloaded tarball (zip) file and execute `runBioDeg.cmd` in Windows or `runBioDeg.sh` in Linux. By doing this you see BiODEG interface showing up on the screen.

## 3.2 Advanced installation for improved performance/flexibility

Building BiODEG and required libraries from source code will increase the performance since the program and the libraries will get optimized for the platform in which they are going to run. Moreover, this enables users to customize the software in the way they want. Additionally, this is an inevitable aspect if you are going to use BiODEG for development purposes or you want to contribute to it.

### 3.2.1 Compiling and installing external libraries

#### 3.2.1.1 PETSc and Qt

BiODEG uses [PETSc](#) for parallel computing. You may choose to build a customized version of PETSc or use the version that comes with FreeFEM. The version that is bundled with FreeFEM has the following build configuration:

```
--with-debugging=0 COPTFLAGS="-O3 -mtune=native" CXXOPTFLAGS="-O3 -mtune=native"
FOPTFLAGS="-O3 -mtune=native" --with-cxx-dialect=C++11 --with-ssl=0 --with-x=0
--with-fortran-bindings=0 --with-cc=/usr/ --with-scalar-type=complex
--with-blaslapack-include= --with-blaslapack-lib="-llapack -lblas"
--with-scalapack --with-metis --with-ptscotch --with-suitesparse --with-suitesparse-lib=
"-Wl, -lumfpack -lklu -lcholmod -lbtf -lccolamd -lcolamd -lcamd -lamd -lsuitesparseconfig"
--with-mumps --with-parmetis --with-tetgen --download-slepc --download-hpddm PETSC_ARCH=fc
```

You may need to build your own version if this configuration is not suitable for you. You can find the instruction for building custom version of PETSc [here](#).

BiODEG UI is developed using [Qt](#) so it should be installed on your system if you want to compile BiODEG UI. You can find the installation instruction for various platforms [here](#).

#### 3.2.1.2 FreeFEM

The full build documentation of FreeFEM is available [here](#), but the following steps is what you need to do to build it on any platform. By default, FreeFEM downloads and builds PETSc during the build process.

1. Install required prerequisites

```
$ sudo apt-get install cpp freeglut3-dev g++ gcc gfortran m4 make patch pkg-config
wget python unzip liblapack-dev libhdf5-dev libgsl-dev autoconf automake
autotools-dev bison flex gdb git cmake
$ sudo apt-get install mpich
```

2. Make a new directory for FreeFEM and navigate to it

```
$ cd
$ mkdir FreeFEM
$ cd FreeFEM/
```

3. Clone the source code repository and navigate to the downloaded directory

```
$ git clone https://github.com/FreeFem/FreeFem-sources.git
$ cd FreeFem-sources/
```

4. Generate the configure scripts

```
$ autoreconf -i
```

5. Run the configure script to specify build options, including the location to install the program

```
$ ./configure --enable-download --enable-optim
--prefix=/home/<your_profile>/FreeFEM/freefem-install
```

6. Download the source code of 3rd-party libraries

```
$ ./3rdparty/getall -a
```

7. Build PETSc and all the 3rd-party libraries

```
$ cd 3rdparty/ff-petsc/
$ make petsc-slepc
```

8. Navigate back and reconfigure the build

```
$ cd -
$ ./reconfigure
```

9. Build the source code of FreeFEM using 4 parallel processes (or nay other number you like

```
$ make -j4
```

10. Check the build by running some examples

```
$ make -j2 check
```

11. Install the built binaries to the specified directory

```
$ make install
```

12. Navigate to the installation location and run FreeFEM

```
$ cd ../freefem-install/  
$ cd bin/  
$ ./FreeFem++
```

13. Navigate to the home directory and add FreeFEM to the PATH variable in the `.bashrc` file

```
$ cd  
$ nano .bashrc
```

and add `export PATH=$PATH:/home/<your_profile>/FreeFEM/freefem-install/bin` to the end of the file and save (press *Ctrl+X* and then *Y*).

After doing this, you should be able to run FreeFEM. Start a new terminal and run `FreeFem++` and `FreeFem++-mpi`. Seeing no error in the output means that you have successfully installed it.

### 3.2.2 Building and installing BioDeg

Since BIODEG UI is developed using Qt, compiling the source files is quite straightforward. Upon installing Qt on your machine, clone the [BIODEG UI repository](#) and follow one of the following scenarios to build it.

#### 3.2.2.1 Build BioDeg UI using Qt Creator IDE

This is the simplest technique to build the program, and it has a similar procedure for all the supported platforms. Qt Creator is the default IDE for Qt development, so it is automatically installed along with Qt. Simply open the Qt project file (`CMakeLists.txt`) in Qt Creator (by executing `qtcreator CMakeLists.txt` or selecting *File->Open Project* from the IDE) and build the project (by pressing *Shift+B*).

#### 3.2.2.2 Build BioDeg UI using Qt tools

Building the source files using CMake is also quite simple. Navigate to the source files directory (the cloned repository) and run the following commands (this assumes that you have already added Qt `bin` directory to the `'PATH'` variable so that the CMake script can find Qt libraries and binaries):

```
$ mkdir build  
$ cd build  
$ cmake ..  
$ make
```

In Windows, you may need to call the correct build system installed along with Qt. For example, by assuming that you have installed the MinGW integration for Qt, the `make` command should be written as `mingw32-make`. Moreover, in this case, you need to call CMake with a suitable generator, so the `cmake` command should be replaced by something like `cmake -G "MinGW Makefiles"` (don't forget to insert the double dots).

After doing this, you can find the BIODEG UI executable in the source directory and run it by executing `./BioDeg` in Linux or `.\BioDeg.exe` in Windows.

## 4 Running BioDeg

### 4.1 Configuring the simulation

#### 4.1.1 Example 1

#### 4.1.2 Example 2

## 5 Future extensions to BioDeg

The future versions of BioDeg will focus on implementing the following methodologies/features.

- Extending the core models to capture more complex chemistry of biodegradation by considering more reactions occurring in buffered solutions.
- Adding support for more base materials like Fe and Zn.
- Considering the effect of alloying elements and complex compositions.
- Adding more post-processing features to BIODEG UI.
- Adding basic visualization to BIODEG UI using ParaView Glance.
- Improving the performance of fluid flow solver by employing a gradient-based solver.
- Considering GPU support by enabling GPU computing in recent versions of PETSc.

## 6 Finding answers to more questions

If you have questions that go beyond this manual, there are a number of resources:

- For questions/suggestions about BIODEG installation, bugs, or similar stuff please use the [BIODEG issue tracker](#).
- BIODEG is primarily based on the [FreeFEM](#). If you have particular questions about FreeFEM, contact the community at <https://community.freefem.org/>.
- If you have specific questions about BIODEG that are not suitable for public and archived mailing lists, you can contact the primary developer and mentor:
  - Mojtaba Barzegari: [mojtaba.barzegari@kuleuven.be](mailto:mojtaba.barzegari@kuleuven.be).
  - Liesbet Geris: [liesbet.geris@kuleuven.be](mailto:liesbet.geris@kuleuven.be) (Mentor).

## A Run-time input parameters

The underlying description of the input parameters also includes a “Standard/Advanced” label, which signifies whether an input parameter is a standard one or an advanced level parameter. The default values of the “Advanced” parameters are good enough for almost all cases. However, in some cases user may need to



use “Advanced” labeled parameters. For user convenience, all input parameters are also indexed at the end of this manual in Section [Index](#).

## A.1 Geometry and mesh parameters

- *Parameter name:* `import_mesh`

*Default:* true (1)

*Description:* [Standard] Boolean parameter specifying whether an external mesh file is imported or a container box as well as a cubic scaffold would be created on the fly for simulation.

*Possible values:* A boolean value (1 or 0)

- *Parameter name:* `mesh_file`

*Default:* Should be provided

*Description:* [Standard] Path to the input mesh file (in MEDIT `.mesh` format), which can be either absolute or relative. Is relevant only if parameter `import_mesh` is set to TRUE.

*Possible values:* Any string value

- *Parameter name:* `label_scaffold`

*Default:* 1

*Description:* [Standard] The label of the (volume) region supposed to be scaffold in the input mesh (can be viewed in programs like GMSH before importing into BIODEG).

*Possible values:* Any positive integer value

- *Parameter name:* `label_medium`

*Default:* 2

*Description:* [Standard] The label of the (volume) region supposed to be the medium (electrolyte) in the input mesh.

*Possible values:* Any positive integer value

- *Parameter name:* `label_wall`

*Default:* 3

*Description:* [Advanced] The label of the surface in the input mesh to be assigned as wall (no slip boundary condition) in the fluid flow simulations.

*Possible values:* Any positive integer value

- *Parameter name:* `label_inlet`

*Default:* 4

*Description:* [Advanced] The label of the surface in the input mesh to be assigned as flow inlet (constant velocity boundary condition) in the fluid flow simulations.

*Possible values:* Any positive integer value

- *Parameter name:* `label_outlet`

*Default:* 5

*Description:* [Advanced] The label of the surface in the input mesh to be assigned as flow outlet (zero pressure boundary condition) in the fluid flow simulations.

*Possible values:* Any positive integer value

- *Parameter name:* `box_length`

*Default:* 20.0

*Description:* [Standard] In case of `import_mesh` being FALSE, specifies the length of the container box (for the electrolyte) in mm.

*Possible values:* Any positive floating point number

- *Parameter name:* `cube_size_x`

*Default:* 13.0

*Description:* [Standard] In case of `import_mesh` being FALSE, specifies the length of the scaffold cuboid along the x axis in mm.

*Possible values:* Any positive floating point number

- *Parameter name:* `cube_size_y`

*Default:* 13.0

*Description:* [Standard] In case of `import_mesh` being FALSE, specifies the length of the scaffold cuboid along the y axis in mm.

*Possible values:* Any positive floating point number

- *Parameter name:* `cube_size_z`

*Default:* 4.0

*Description:* [Standard] In case of `import_mesh` being FALSE, specifies the length of the scaffold cuboid along the z axis in mm.

*Possible values:* Any positive floating point number

- *Parameter name:* `mesh_size`

*Default:* 32

*Description:* [Standard] Number of elements on each edge of the container box, so a higher number means a finer mesh. The mesh size of the cuboid will be adjusted accordingly or can be adaptively refined by setting parameter `refine_initial_mesh` to TRUE.

*Possible values:* Any positive integer number

- *Parameter name:* `refine_initial_mesh`

*Default:* false (0)

*Description:* [Advanced] A boolean parameter specifying if the mesh (no matter if imported or generated) should be adaptively refined on the metal-medium interface (corrosion surface). This affects the beginning of the simulation only (on the initial mesh).

*Possible values:* A boolean value (1 or 0)

- *Parameter name:* `mshmet_error`

*Default:* 0.01

*Description:* [Advanced] Since the open source tool `mshmet` is used for creating a metric for refining the mesh on the level set signed distance function, a tolerance should be specified for it. A lower value results to a finer mesh.

*Possible values:* Any floating point number

- *Parameter name:* `mesh_size_min`

*Default:* 0.04

*Description:* [Advanced] Specifies the smallest element size to be passed to the `tetgen` mesh generator for refining the initial mesh.

*Possible values:* Any floating point number

- *Parameter name:* `mesh_size_max`

*Default:* 0.8

*Description:* [Advanced] Specifies the largest element size to be passed to the `tetgen` mesh generator for refining the initial mesh.

*Possible values:* Any floating point number

## A.2 Materials and boundary conditions parameters

- *Parameter name:* `material_density`

*Default:* 1.735e-3

*Description:* [Standard] The density of the metallic material. The default value is the density of magnesium.

*Possible values:* Any floating point number

- *Parameter name:* `film_density`

*Default:* 2.3446e-3

*Description:* [Standard] The density of the protective film that forms on the corrosion surface. The default value is the density of magnesium hydroxide.

*Possible values:* Any floating point number

- *Parameter name:* `material_satur`

*Default:* 0.134e-3

*Description:* [Advanced] The saturation concentration at which the metallic material (here, the ions) saturates through the medium. The default value is defined for magnesium ions.

*Possible values:* Any floating point number

- *Parameter name:* `material_eps`

*Default:* 0.55

*Description:* [Advanced] The porosity of the formed protective layer in the range [0, 1].

*Possible values:* Any floating point number between 0 and 1

- *Parameter name:* `material_tau`

*Default:* 1.0

*Description:* [Advanced] The tortuosity of the formed protective layer.

*Possible values:* Any floating point number

- *Parameter name:* `d_mg`

*Default:* 0.05

*Description:* [Standard]

*Possible values:* The diffusion coefficient of the metallic ions transport. This parameter is one of the most effective ones on the rate of degradation.

- *Parameter name:* `d_cl`

*Default:* 0.05

*Description:* [Standard] The diffusion coefficient of the chloride ions transport

*Possible values:* Any floating point number

- *Parameter name:* `d_oh`

*Default:* 25.2

*Description:* [Standard] The diffusion coefficient of the hydroxide ions transport

*Possible values:* Any floating point number

- *Parameter name:* `k1`

*Default:* 7.0

*Description:* [Standard] The reaction rate at which the chemical reaction of the protective film formation occurs.

*Possible values:* Any floating point number

- *Parameter name:* `k2`

*Default:* 1e15

*Description:* [Standard] The reaction rate at which the chemical reaction of the protective film dissolution occurs.

*Possible values:* Any floating point number
- *Parameter name:* `fluid_nu`

*Default:* 0.85

*Description:* [Advanced] The effective viscosity of the fluid used to simulate hydrodynamics condition.

*Possible values:* Any floating point number
- *Parameter name:* `fluid_in_x`

*Default:* 0.1

*Description:* [Advanced] The X component of the fluid velocity defined on inlet (see `label_inlet`) as the boundary condition for the fluid flow.

*Possible values:* Any floating point number
- *Parameter name:* `fluid_in_y`

*Default:* 0

*Description:* [Advanced] The Y component of the fluid velocity defined on inlet (see `label_inlet`) as the boundary condition for the fluid flow.

*Possible values:* Any floating point number
- *Parameter name:* `fluid_in_z`

*Default:* 0

*Description:* [Advanced] The Z component of the fluid velocity defined on inlet (see `label_inlet`) as the boundary condition for the fluid flow.

*Possible values:* Any floating point number
- *Parameter name:* `initial_cl`

*Default:* 5.175e-6

*Description:* [Standard] Initial concentration of chloride ions in the medium.

*Possible values:* Any floating point number
- *Parameter name:* `initial_oh`

*Default:* 1e-7

*Description:* [Standard] Initial concentration of hydroxide ions in the medium, used for computing pH. The default value indicates a pH of 7.

*Possible values:* Any floating point number

### A.3 Solver parameters

- *Parameter name:* `solve_mg`

*Default:* true (1)

*Description:* [Standard] Boolean parameter indicating whether the equation for material dissolution and ions release should be solved or not. This equation is the most essential equation and in most use-cases should be solved.

*Possible values:* Any boolean value (1 or 0)

- *Parameter name:* `solve_film`

*Default:* true (1)

*Description:* [Standard] Boolean parameter indicating whether the equation for the protective film formation should be solved or not.

*Possible values:* Any boolean value (1 or 0)

- *Parameter name:* `solve_cl`

*Default:* true (1)

*Description:* [Standard] Boolean parameter indicating whether the equation for the transport of chloride ions should be solved or not.

*Possible values:* Any boolean value (1 or 0)

- *Parameter name:* `solve_oh`

*Default:* true (1)

*Description:* [Standard] Boolean parameter indicating whether the equation for the transport of hydroxide ions should be solved or not. This equation is essential for calculating the pH changes, if desired.

*Possible values:* Any boolean value (1 or 0)

- *Parameter name:* `solve_ls`

*Default:* true (1)

*Description:* [Standard] Boolean parameter indicating whether the level set surface tracking equation should be solved or not. Surface tracking is essential for computing the mass loss.

*Possible values:* Any boolean value (1 or 0)

- *Parameter name:* `solve_fluid`

*Default:* false (0)

*Description:* [Standard] Boolean parameter indicating whether the fluid flow equation should be solved or not.

*Possible values:* Any boolean value (1 or 0)

- *Parameter name:* `solve_full_ns`

*Default:* true (1)

*Description:* [Advanced] Boolean parameter specifying which fluid flow equation to solve: the full transient Navier-Stokes equations or a steady-state Stokes equation. A true value (1) results in BIODEG solving the former equation.

*Possible values:* Any boolean value (1 or 0)
- *Parameter name:* `write_fluid_output`

*Default:* true (1)

*Description:* [Advanced] Boolean parameter specifying whether the fluid flow quantities (velocity field components and pressure) should be saved in the output VTK file or not. Requires `write_vtk` be TRUE.

*Possible values:* Any boolean value (1 or 0)
- *Parameter name:* `solve_fluid_each`

*Default:* 10

*Description:* [Advanced] Determines the number of time steps to skip before solving the specified fluid flow equation. For example, if the value is set to 10 (default value), the fluid equation gets solved in time steps 1, 11, 21, ....

*Possible values:* Any positive integer number
- *Parameter name:* `time_step`

*Default:* 0.025

*Description:* [Advanced] The time step value for the numerical schemes of the simulations.

*Possible values:* Any floating point number
- *Parameter name:* `final_time`

*Default:* 21.0

*Description:* [Standard]

*Possible values:* Any floating point number
- *Parameter name:* `do_redistance`

*Default:* true (1)

*Description:* [Advanced]

*Possible values:* Any boolean value (1 or 0)
- *Parameter name:* `redistance_time`

*Default:* 1.0

*Description:* [Advanced]

*Possible values:* Any floating point number

## A.4 Output parameters

- *Parameter name:* `text_output_file`

*Default:* "output/result.txt"

*Description:* [Standard] Path to the text file in which text output, like the mass loss and evolved hydrogen production, are written over time. The path can be relative or absolute.

*Possible values:* Any string value referring to a valid path

- *Parameter name:* `write_vtk`

*Default:* true (1)

*Description:* [Standard] Indicates whether VTK output (in the VTU format) should be written or not. This is required for further post-processing of the results using ParaView.

*Possible values:* Any boolean value (1 or 0)

- *Parameter name:* `vtk_output_name`

*Default:* "output/output"

*Description:* [Standard] The naming scheme for the VTK output. This is mainly for the PVD file, and the final name of the rest of the files will be determined by the number of employed MPI computing nodes and the time step. The number of VTU files saved per time step equals to the number of employed MPI nodes.

*Possible values:* Any string value referring to a valid path

- *Parameter name:* `save_each`

*Default:* 0.25

*Description:* [Standard] The interval of saving results, text and VTK (if selected), to disk.

*Possible values:* Any floating point number

- *Parameter name:* `save_last_state`

*Default:* true (1)

*Description:* [Standard] Indicates whether the last state of the system should be saved or not. The last state will be always saved on a global (non-partitioned) mesh, meaning that it will be a single VTU file, in contrast to normal save in which the number of written files equals to the number of computing nodes.

*Possible values:* Any boolean value (1 or 0)

- *Parameter name:* `output_per_area`

*Default:* false (0)

*Description:* [Advanced] Indicated whether the side hydrogen evolution results should be computed per unit area of exposed surface.

*Possible values:* Any boolean value (1 or 0)



- *Parameter name:* `save_multiplier`

*Default:* 1.0

*Description:* [Advanced] In case of symmetrical conditions, this parameter can be used to multiply the obtained quantitative results to have easier comparison with experimental results. A value of 1.0 indicates no multiplication.

*Possible values:* Any floating point number

- *Parameter name:* `export_scaffold`

*Default:* false (0)

*Description:* [Advanced] Indicates if the degrading material should be exported as a single entity for further analysis like in a structural mechanics simulation. It works based on Mmg level set meshing.

*Possible values:* Any boolean value (1 or 0)

- *Parameter name:* `export_scaffold_each`

*Default:* 1.0

*Description:* [Advanced] In case `export_scaffold` being true, this parameter determines the interval of saving the material mesh to the disk.

*Possible values:* Any floating point number

- *Parameter name:* `export_scaffold_volume`

*Default:* true (1)

*Description:* [Advanced] In case `export_scaffold` being true, this parameter indicates whether a volume mesh should be saved as output or not.

*Possible values:* Any boolean value (1 or 0)

- *Parameter name:* `export_scaffold_surface`

*Default:* true (1)

*Description:* [Advanced] In case `export_scaffold` being true, this parameter indicates whether a surface mesh should be saved as output or not.

*Possible values:* Any boolean value (1 or 0)

- *Parameter name:* `save_initial_mesh`

*Default:* false (0)

*Description:* [Advanced] Determines whether the initial computational mesh should be saved for further debugging or not. Can be used to investigate the mesh refinement if it is asked by setting `refine_initial_mesh` to true.

*Possible values:* Any boolean value (1 or 0)

- *Parameter name:* `save_initial_partitioned_mesh`

*Default:* false (0)

*Description:* [Advanced] Determines whether the partitioned computational mesh should be saved for further debugging or not. Can be used to debug and view the output of the mesh partitioning process before going for the actual simulation.

*Possible values:* Any boolean value (1 or 0)

## Index of run-time parameters with section names

The following is a listing of all run-time parameters, sorted by the section in which they appear.

### Geometry and mesh

- box length, [10](#)
- cube size x, [10](#)
- cube size y, [10](#)
- cube size z, [10](#)
- import mesh, [9](#)
- label inlet, [9](#)
- label medium, [9](#)
- label outlet, [10](#)
- label scaffold, [9](#)
- label wall, [9](#)
- mesh file, [9](#)
- mesh size, [10](#)
- mesh size max, [11](#)
- mesh size min, [11](#)
- mshmet error, [11](#)
- refine initial mesh, [10](#)

### Materials and boundary conditions

- d cl, [12](#)
- d mg, [12](#)
- d oh, [12](#)
- film density, [11](#)
- fluid in x, [13](#)
- fluid in y, [13](#)
- fluid in z, [13](#)
- fluid nu, [13](#)
- initial cl, [13](#)
- initial oh, [13](#)
- k1, [12](#)
- k2, [13](#)
- material density, [11](#)

- material eps, [12](#)

- material satur, [12](#)

- material tau, [12](#)

### Output

- export scaffold, [17](#)
- export scaffold each, [17](#)
- export scaffold surface, [17](#)
- export scaffold volume, [17](#)
- output per area, [16](#)
- save each, [16](#)
- save initial mesh, [17](#)
- save initial partitioned mesh, [18](#)
- save last state, [16](#)
- save multiplier, [17](#)
- text output file, [16](#)
- vtk output name, [16](#)
- write vtk, [16](#)

### Solver

- do redistance, [15](#)
- final time, [15](#)
- redistance time, [15](#)
- solve cl, [14](#)
- solve film, [14](#)
- solve fluid, [14](#)
- solve fluid each, [15](#)
- solve full ns, [15](#)
- solve ls, [14](#)
- solve mg, [14](#)
- solve oh, [14](#)
- time step, [15](#)
- write fluid output, [15](#)